

Non-invasive, drug free  
therapy for Sports



# SCENAR *for Sports*

[www.scenar.com.au](http://www.scenar.com.au)

Self Controlled Energo Neuro Adaptive Regulator

Designed for Professional athletes

Sports people

Coaches

Fitness instructors

Home users



# Try SCENAR™ SPORT - Reduce Trauma Downtime

**RITMSCENAR™ Sport** is an electro-therapy device with built-in electrode that delivers complex electrical impulses through the skin. These impulses stimulate the Nervous System to generate Neuropeptides, the key biochemicals needed by the body to heal itself.

**Neuropeptides may be considered as our body's own internal pharmacy.**

**RITMSCENAR™ Sport** controls and automatically changes the impulse waveforms and frequencies to prevent the body from adapting to the stimulation.

**RITMSCENAR™ Sport** is designed for sportsmen to **increase performance effectiveness, reduce trauma downtime** and **speed-up rehabilitation process.**

In Australia, the SCENAR devices are intended to be used as an analgesic peripheral nerve stimulator, to treat pain associated with surgery, trauma, musculoskeletal problems, bursitis or dental problems. It can also have further applications in facilitating functional restoration and physical therapy.

**RITMSCENAR™ Sport** is a **portable, light, easy to operate** device that can be used at competitions, training sessions, travelling or other locations of a sports team.

**Two model variants are available:**

RITMSCENAR™ **Sport** with Light Emitting diodes  
RITMSCENAR™ **Sport D** with LCD display.





"SCENAR™ supported my nervous system, encouraged my body to heal, allowed me to maintain strength and fitness so I could compete at my optimum level "

**Layne Beachley**, Australia  
7X ASP World Surfing Champion

"I added SCENAR™ treatment to my recovery plan. Without SCENAR™ treatment I wouldn't have been able to play at the 2009 NRL Grand Final"

**Nathan Cayless**, Australia  
Sydney NRL Parramatta Eels Captain



"SCENAR™ therapy benefits appear to be almost limitless for health and wellbeing."

**Torah Bright**, Australia  
Olympic Gold medalist  
TTR World Snowboard Champion

## Device Parameters

### RITMSCENAR™ Sport / Sport D

Working Modes	Subjective (Manual) Dosed (Automatic)
Dose indication	Sound and light/Symbolic
Amplitude modulation	3 : 1
Frequency modulation	30 -120 Hz, cycle 7 sec
Influencing impulse frequencies	14, 60, 90, 340 Hz
Intensity <i>(the number of impulses in a packet)</i>	1
Amplitude of the stimulating pulse at a standard load	(1,7-2,8) ... (100-140) V
Indicator	LED / LCD
Power supply	one 9V alkaline battery
Temperature	10 to 35C (Operating) -50 to 50C (Storage & Transportation)
Weight	0.2kg
Dimensions	130 x 50 x 25mm
Warranty	2 years

Certificates:



Australia  
ARTG  
ID: 140659

U.S.A.  
FDA  
K092117

[www.scenar.com.au](http://www.scenar.com.au)

**Branch office for Asia Pacific region**

**RITM Australia Pty Ltd**

7 Belsham Road  
Kariong NSW 2250  
Australia

Tel: +61 2 8011 4217  
Fax: +61 2 9423 4681  
info@scenar.com.au  
www.scenar.com.au

**TRAINING**

Additional training for RITMSCENAR™ Sport users is available. It covers theoretical and practical applications of the device and technology.

For more information about training options, please contact the Authorized Training Centre:

**RITM SCENAR Institute Pty Ltd.**

Tel: +61 2 8005 6461  
info@scenarinstitute.org  
www.scenarinstitute.org

Distributor:

Self Controlled Energo Neuro Adaptive Regulator

**SCENAR**<sup>TM</sup> **Sport**