

RITM

SCENAR[®]

PERSONAL DEVICES

www.scenar.com.au



*Non-invasive, drug free
therapy for pain relief*

RITM
SCENAR AUSTRALIA



RITMSCENAR™ Sport D, Sport & Home


Sport D, Sport and Home, are electro-therapy devices with built-in electrodes that deliver complex electrical impulses through the skin to give pain relief.

SCENAR™ interaction

The SCENAR™ controls and automatically changes the impulse waveforms and frequencies to prevent the body from adapting to the stimulation.

RITMSCENAR™ usage

RITMSCENAR personal range of devices are designed for home users to manage their pain and to speed up the recovery process from acute and chronic pain, sport injury and post operative pain. The personal devices are light, compact and can be used almost anywhere.



"With consistent SCENAR treatment and maintenance, I was able to avoid surgery and retirement and return to competitive surfing without pain"

Layne Beachley, Australia
7x ASP Women's World
Surfing Champion




"Since 2009 SCENAR therapy is part of my daily injury prevention and management."

Nathan Cayless, Australia
Former NRL Parramatta
Eels Captain



"I am using RITM SCENAR for both injury and well-being while competing and travelling"

Torah Bright, Australia
World Snowboard Champion
and Olympic Gold medallist

FUNCTION/DEVICE	RITMSCENAR Sport D	RITMSCENAR Sport	RITMSCENAR Home
			
Indication	LCD display	LED	LED
Skin contact indication	Yes	No	No
Dose Modes	Dose 1 and 2	Dose 1 and 2	Dose 1 and 2
Amplitude modulation	“3 : 1”	“3 : 1”	“3 : 1”
Frequency modulation	30Hz – 120Hz	30Hz – 120Hz	30Hz – 120Hz
Influencing frequencies	14Hz, 60Hz, 90Hz, 340Hz	14Hz, 60Hz, 90Hz, 340Hz	90Hz
Power supply	9V battery	9V battery	9V battery
Weight	0.2kg	0.2kg	0.2kg
Dimensions	13 x 5 x 2,5cm	13 x 5 x 2,5cm	13 x 5 x 2,5cm
Warranty	2 years	2 years	2 years


BRANCH OFFICE FOR ASIA PACIFIC REGION
RITM AUSTRALIA PTY LTD

P O Box 615 St Ives NSW 2075 Australia
Tel: +61 2 8011 4217 • Fax: +61 2 9423 4681
info@scenar.com.au • www.scenar.com.au
training web site • www.ritm-scenar.com.au

DISTRIBUTOR



TRAINING



All RITMSCENAR® personal devices come with a comprehensive User's Manual and free one-year online course. The training resources include information about how to operate with the device, principles of SCENAR technology, practical treatment protocols of most common pain complaints.

The online training provides a dedicated forum for questions and discussions about SCENAR devices and treatment techniques.



Our online training web site is:
www.ritm-scenar.com.au.

Additionally, we organise from time to time live seminars for Home Users. Basic and Advanced levels are available in Australia and New Zealand.